
ELECTRO DETACH GENERATOR V4

Detach Point Detection System

The Electro Detach Generator v4 detects and indicates by sound and light signals the best Detach Point even when the position of the second marker of a microcatheter is not detectable under a X-ray fluoroscopy, Easy to use, disposable, sterilized, detachment of the coil can be done in single step operation by pressing 2 buttons simultaneously.



stays within the microcatheter.

The green lamp is lit when the detaching part of the ED COIL has come out of the microcatheter.

Orange lamp The orange lamp is lit when a current leak or other trouble







Generator v4 is powered on.

DETACH button Detaching output is generated when the left and right buttons are pressed simultaneously while the

Manual Detachment -

When it is necessary to detach the ED COIL while the red lamp is lit, press and hold the left and right buttons simultaneously for longer than 2 seconds, and a detaching output is generated.

If the green lamp is lit while the detaching part is still within the microcatheter, press the reset button to restart detecting the detach point.



The red lamp is lit when the ED clip is connected to the pusher.




The green lamp is lit when the detaching part of the ED COIL has come out of the microcatheter. Subsequently, detachment of the coil can be done by pressing 2 buttons simultaneously





Patient-side clip to be used with Hypodermic Needle made of stainless-steel without







Do not re-use STERILEEO Sterilized using ethylene oxide



Consult instructions for use

KANEKA CORPORATION CONTACT TEL. +81-3-5574-8136

TEL. +81-6-6226-5256

KANEKA PHARMA EUROPE N.V.







ED COIL—Electro Detach Coil—

"Reliable Coil Embolization" achieved by "Accurate Coil Placement and Tight Packing"





Instant Detachment


A platinum coil is connected to the pusher with a short poly vinyl alcohol (PVA) rod. The PVA rod is thermally melt down (at around 70°C) by high-frequency current generated by the Electro Detach Generator v4, and the coil is instantly detached.

Detach Point Detection

The Electro Detach Generator v4 detects and indicates by sound and light signals the best Detach Point even when the position of the second marker of a microcatheter is not detectable under a X-ray

Soft Tip Pusher


The tip (through 30 mm from the distal end) is very flexible, therefore, it prevents kick-back of the microcatheter while deploying the coil in the aneurysm



ED COIL10 ExtraSoft


Soft and flexible coil for finishing

Extremely flexible coil that facilitates tight packing





2-4ES (φ2mm x 4cm) tightly packed in a 2 mm aneurysmal model.



anti-stretching characteristics.

The stretch resistant wire provides flexible and


MOMMAN

EIOCOIL 10 [Infini]

A unique large diameter coil provides flexible filling regardless of the size and shape of an aneurysm.

- · Conforming outward to reduce compartment space
- · Compactly folding and entering into even a gap space
- · Conforms in accordance with the shape of an aneurysm even with irregular shape
- · Provides a tight embolization with a small number of coils for a vessel occlusion

Waving Inner Wire Technology



Waving Inner Wire Technology enables coil-folding in smaller size despite its large



Infini Soft's large coil-diameter reduces the risk of coil-protrusion through a stent strut into the parent artery, and it enables its use for both framing and filling



aneurysm wall, thereby contributing to uniform packing even in an aneurysm


ED COIL14 Standard Office

Helical coil spreading randomly

- The coil with 0.014 inch primary coil-diameter may provide sufficient coil embolization with a small number of coils
- The unique aSpiral coil structure provides random coil movement by deflection according to the aneurysm wall
- · Strong framing by the unique aSpiral structure
- · Suitable for both framing and filling
- ≤ 6 mm (small pitch): suitable for filling
- ≥ 7 mm (large pitch): suitable for both framing and filling



The unique aSpiral shape, designed to form a strong frame with random loops by having 1:1 ratio between coil diameter and pitch (≥ 7 mm coil)



The random coil spreading with the αSpiral structure enables to change direction of the coil expanding in accordance with the shape of the wall of



αSpiral (φ10mm×L30cm) in a model of a 10 mm aneurysm